



Superior Scleral Lens Design
Gives You the Ability to Fit
Irregular and Regular Corneas



IRREGULAR CORNEAS

16.3 mm

Flexibility

- FlexFit's Auto-Adjusting design allows you to order custom diameters from 15.5mm to 17.0mm in .1mm increments while maintaining desired apical clearance
- Incremental adjustments to the fit can be accomplished through changes to the Central Clearance, Peripheral Corneal Clearance, Limbal Clearance and Scleral Landing Zones
- Available in a front surface toric design

Dual Depth Periphery™

This feature allows ICD FlexFit™ to land appropriately on the sclera, delivering:

- Rotational stability
- Extended wear time
- Improved lens comfort

REGULAR CORNEAS

14.8 mm

Auto Adjusting Features

- FlexFit diameter from 14.5mm to 15.5mm in .1mm increments to maintain apical clearance
- Innovative scribe mark technology on trial lenses ensures quick, accurate limbal clearance



ICDLens.com