



eaglet-eye

appenzeller
kontaktlinsen

Eye Surface Profiler and i-Shape


A MOST POWERFUL COMBINATION

Did you know that only 6% of scleras are spherical and two-thirds are asymmetrical or irregularly shaped? Biometric lenses are custom-tailored to the exact geometry of the eye and offer the highest comfort and vision, even in the most complex cases.



CORNEO-SCLERAL PROFILOMETRY

- **Save time** with quick measurements at the push of a button.
- **No more guesswork** – real data without extrapolation.
- **Extensive eye coverage** with a scanning range of up to 20 mm.
- **Precision** down to the detail with over 500,000 captured data points.
- Ideal for both regular and **irregular eyes**.



PERSONALIZED CONTACT LENSES

- **Software-based** fitting without trial lenses.
- **Improved fitting accuracy.**
- Manufacture of a **freeform scleral lens** based on individual eye profiles using ESP data.
- **More comfortable** thanks to optimal landing and pressure distribution on the sclera.
- **Reduced** time and cost effort.
- **Fits any eye**, no matter how irregular.



Included with your Eye Surface Profiler (ESP):

- First Lens Fit software for all lens modalities
- 2 years of free software updates
- Online training levels 1 & 2
- Lifetime access to the Eaglet Eye Academy
- Comprehensive optometry support

Your benefits with i-Shape®:

- 50% discount on an i-MATRIX fitting set
- Voucher for 1 pair of i-Shape® scleral lenses
- Seminar bundle online course FIT1 to FIT5
- Consultation with your Appenzeller contact lens experts

* Exclusive offer for the purchase of an Eye Surface Profiler between November 1, 2024, and January 31, 2025 (including Opti). Subject to change.